

# Canespro Fungal Nail Treatment Set



## Category

Footcare

## Manufacturer

BAYER PLC

## Pack size

10g

## RRP

£38.29

## Legal Status

MD



## Uses

Removes the infected parts of the nail in 2-3 weeks.



## Format and Ingredients

Treatment set containing:

10g tube Canespro Urea 40% ointment, 22 waterproof plasters, 1 plastic nail scraper,



## Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.

Step 2: Cover the treated nail with one of the plasters provided and keep covered for 24 hours.

Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.

Repeat steps 1, 2 and 3 daily for 2-3 weeks.