


# Canespro Fungal Nail Treatment Set



**Category**  
Footcare  
**Manufacturer**  
BAYER PLC  
**Pack size**  
10g  
**RRP**  
£38.29  
**Legal Status**  
MD



## Uses

Removes the infected parts of the nail in 2-3 weeks.



## Format and Ingredients

Treatment set containing:  
10g tube Canespro Urea 40% ointment, 22 waterproof plasters, 1 plastic nail scraper,



## Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.  
Step 2: Cover the treated nail with one of the plasters provided and keep covered for 24 hours.  
Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.  
Repeat steps 1, 2 and 3 daily for 2-3 weeks.